


REVISIONS		
DATE	REVISOR	DESCRIPTION
		Rev



Laser Engraved
Serial # and fiber core size

Material: Arcap (Stainless Steel)
Finished Ferrule length: .3853" to .3865"

DRAWN B. Rogers	DATE 07/23/04	TOLERANCES PRIMARY UNITS ARE INCHES AND SECONDARY UNITS ARE MILLIMETERS UNLESS OTHERWISE SPECIFIED	380 MAIN STREET DUNEDIN, FL 34698 (727) 733-2447 <i>Ocean Optics Inc.</i>			TITLE OSMA data sheet
PROJECT ENGINEER B. Rogers		INCHES DECIMAL: 0.005 FRACTIONS: 1/32 ANGLE: 4:3	MILLIMETERS 0.25 \pm 0.10 ANGLE: 30:5	SCALE NTS	SHEET NO 1 of 1	DWG. NO.
DO NOT SCALE DRAWING			<small>THE MATERIAL AND THE CONTENTS HEREOF IS CONSIDERED PROPERTY AND CONFIDENTIAL INFORMATION OF OCEAN OPTICS, INC. DISCLOSURE TO UNAUTHORIZED PERSONS OR ORGANIZATIONS, REPRODUCTION, OR OTHER USE IS PROHIBITED WITHOUT PRIOR WRITTEN CONSENT BY OCEAN OPTICS, INC.</small>			Rev A